



## Escape-Carry Chair ST-PLUS



### Evacuation ascending the stairs, while in a seated position.

#### What's an Escape-Carry Chair?

An **Escape-Carry Chair** is the solution to evacuate people to safety in a seated position up and/or down a flight of stairs. A simple procedure can make the **Escape-Chair** ready within a few seconds.

#### How does an Escape-Carry Chair work?

The unique glide and brake system makes the use of the **Escape-Carry Chair** possible on all different kinds of stairs, from wood to concrete. The braking system is optimized by the weight of the evacuee, in combination with the pressure of the operator pushing down on the handle in the direction of the stairs.

Ascending evacuation should always be operated by at least 2 people. Ensure to make use of the support braces.

#### Advantages of the Escape-Carry Chair ST-PLUS:

- Used for ascending and/or descending evacuation on the stairs.
- Ideal for quick evacuation.
- Descending evacuation can be operated by one person.
- Speed of descent to be determined by the user.
- Ready to use within seconds.
- Can easily manoeuvre into the evacuation flow.
- Forms minimal blockage on the stairs.
- Equipped with carrying handles.
- Fits into your emergency plan.
- Safe closed gliding belt cassette.
- CE, TÜV & ISO 13485 certified.

Specifications	Escape-Carry Chair ST-PLUS
Article number:	10120
Dimensions folded up (height x width x depth):	120 x 51 x 21 cm
Maximum payload:	225 kg
Maximum stair gradient:	40°
Weight:	± 14 kg
Ergonomic backrest:	double
Demonstration and advice provided:	✓
Ready for use within seconds:	✓
Speed of descent to be determined by the user:	✓
Safe closed gliding belt cassette:	✓
Reflective strips:	✓
CE, TÜV & ISO 13485 certified:	✓
Delivery included:	✓
4-language manual:	✓
Identity sign:	✓
Wall hooks:	✓
Footrest:	optional
Armrests:	optional
Padded (upholstered) headrest:	✓
Comfort seat:	✓
Two height adjustable operating handles:	✓
Warranty on non-wearing parts:	10 years

#### Our advice:

Make sure you have an **Escape-Carry Chair** available in your organisation so that you're prepared for any emergency, for up and/or down evacuation.

**"Failure to prepare for an emergency has its consequences."**



Escape Mobility USA.

1261 S Lyon Street Suite 405 Santa Ana, CA 92705, USA

T +1 (855) 470-4410

W [www.escape-mobility-usa.com](http://www.escape-mobility-usa.com)

