



# Escape-Chair ST-Plus



## Evacuation ascending the stairs, while in a seated position.

### What's an Escape-Chair ST-Plus?

An **Escape-Chair ST-Plus** is the solution to evacuate people to safety in a seated position down a flight of stairs. A simple procedure can make the **Escape-Chair ST-Plus** fully functional within a few seconds.

### How does an Escape-Chair work?

The unique glide/brake system makes the braking of the **Escape-Chair ST-Plus** possible on all different kinds of stairs, from wood to concrete. This braking system is optimized by the weight of the evacuee, in combination with the pressure of the operator pushing down on the handle in the direction of the stairs.

### Advantages of the Escape-Chair ST-Plus:

- Ideal for quick evacuation.
- Can be operated by a single person.
- Ready for use within seconds.
- Carrying people is no longer necessary.
- Speed of descent to be determined by the user (extra brake is unnecessary).
- Easily maneuvering into the evacuation flow.
- Forms minimal blockage on the stairs.
- Fits into your emergency plan.
- Safe closed gliding belt cassette.
- CE, TÜV & ISO 13485 certified.

Specifications	Escape-Carry Chair® ST-Plus
Article number:	10040
Dimensions folded up (height x width x depth):	120 x 57 x 21 cm
Maximum payload:	225 kg
Maximum stair gradient:	40°
Weight:	± 13 kg
Ergonomic backrest:	double
Demonstration and advice provided:	✓
Ready for use within seconds:	✓
Speed of descent to be determined by the user:	✓
Safe closed gliding belt cassette:	✓
Reflective strips:	✓
CE, TÜV & ISO 13485 certified:	✓
Delivery included:	✓
4-language manual:	✓
Identity sign:	✓
Wall hooks:	✓
Footrest:	✓
Armrests:	✓
Padded (upholstered) headrest:	✓
Comfort seat:	✓
Two height adjustable operating handles:	✓
Warranty on non-wearing parts:	10 years

### Our advice:

Make sure you have an **Escape-Chair ST-Plus** available in your organisation so that you're prepared for any emergency, for up and/or down evacuation.

**"Failure to prepare for an emergency has its consequences."**



Escape Mobility USA.

1261 S Lyon Street Suite 405 Santa Ana, CA 92705, USA

T +1 (855) 470-4410

W [www.escape-mobility-usa.com](http://www.escape-mobility-usa.com)

