



## Escape-Mattress Compact



### Best choice for healthcare environment and public buildings

#### What's an Escape Mattress?

**Escape-Mattress** models are designed for individuals who would not be able to comfortably get into an evacuation chair, or where evacuation chairs may struggle (for example when evacuating an individual down a spiral staircase). Evacuation mattresses are ideal for medical centres and hospitals, as they cater for complex healthcare needs. E.g. bed-bound and unconscious patients. **The Escape-Mattress** is an "A-label quality" product and produced in-house.

Our mattresses are CE certified & TÜV certified up to 150kg. We guarantee the highest quality standards.

The **Escape-Mattress Compact** is a complete **Escape-Mattress** but its smaller dimensions allows it to be used more flexibly. The **Escape-Mattress Compact** is rolled up and stored in the storage bag that is included. This **Escape-Mattress Compact** is more suitable, so that it can be taken to different situations where an evacuation aid is required.

The **Escape-Mattress** is made of fire resistant material that is resistant to blood, urine, acids and alkaline as well.

CE certified and TÜV tested and certified for 150kg payload.

There is no better and more durable quality available on the market than any of our **Escape-Mattress** range.

#### Advantages of the EscapeMattress Compact:

- Stored double-folded (space saving)
- 4 carrying handles on each side

#### Protection of evacuee/patient:

- 2 body flaps – 1 per side
- Foot bag
- Headrest

#### Specifications Escape Mattress Compact

Article number:	21030
Free demonstration and expert advice in advance:	✓
Ready to use within a few seconds:	✓
Heat and flame resistant:	✓
Dimensions (length x width):	200 x 68 cm
Mattress thickness:	6 cm
Including storage bag or case	✓
Maximum carrying capacity:	150 kg
Space-saving:	Can be rolled up
Weight:	7.5 kg incl. bag
Adjustable headrest with immobilisation band:	✓
Bag dimensions:	80 x 70 x 21 cm
Warranty:	2 years
4-language manual:	✓
TÜV & CE certified	✓

#### Usage:

The **Escape-Mattress Compact** is mainly used in healthcare environments and (public) buildings.



Escape Mobility USA.

1261 S Lyon Street Suite 405 Santa Ana, CA 92705, USA

T +1 (855) 470-4410

W [www.escape-mobility-usa.com](http://www.escape-mobility-usa.com)

