



## Escape-Mattress Fold



### Best choice for healthcare environment and public buildings

#### What's an Escape Mattress?

**Escape-Mattress** models are designed for individuals who would not be able to comfortably get into an evacuation chair, or where evacuation chairs may struggle (for example when evacuating an individual down a spiral staircase). Evacuation mattresses are ideal for medical centres and hospitals, as they cater for complex healthcare needs. E.g. bed-bound and unconscious patients. **The Escape-Mattress** is an "A-label quality" product and produced in-house.

Our mattresses are CE certified & TÜV certified up to 150kg.

We guarantee the highest quality standards.

#### Usage:

The **Escape-Mattress Fold** is mainly used in healthcare environments and (public) buildings.

The **Escape-Mattress Fold** is folded double so that makes it easy to move and store. Once unfolded, this evacuation mattress is immediately ready to use. This mattress can be stored both hygienically and dust-free, by purchasing either a storage bag or an Escape-Locker.

The **Escape-Mattress** is made of fire resistant material that is resistant to blood, urine, acids and alkaline as well.

CE certified and TÜV tested and certified for 150kg payload.

There is no better and more durable quality available on the market than any of our **Escape-Mattress** range.

#### Advantages of the EscapeMattress Fold:

- Stored double-folded (space saving)
- 4 carrying handles on each side

#### Protection of evacuee/patient:

- 4 body flaps – 2 per side
- Foot bag
- Headrest

#### Specifications Escape Mattress Fold

Article number:	21010
Free demonstration and expert advice in advance:	✓
Ready to use within a few seconds:	✓
Heat and flame resistant:	✓
Dimensions (length x width):	200 x 68 cm
Also available in a 55cm wide version:	✓
Mattress thickness:	8 cm
Including storage bag or case:	Optional
Fits in an Escape-Locker:	✓
Maximum carrying capacity:	150 kg
Space-saving:	Double folding
Weight:	7.2 kg
Adjustable headrest with immobilisation band:	✓
Warranty:	2 years
4-language manual:	✓
TÜV & CE Certified:	✓



Escape Mobility USA.

1261 S Lyon Street Suite 405 Santa Ana, CA 92705, USA

T +1 (855) 470-4410

W [www.escape-mobility-usa.com](http://www.escape-mobility-usa.com)

