



## Escape-Mattress Solo



### Best choice for healthcare environment and public buildings

#### What's an Escape Mattress?

**Escape-Mattress** models are designed for individuals who would not be able to comfortably get into an evacuation chair, or where evacuation chairs may struggle (for example when evacuating an individual down a spiral staircase). Evacuation mattresses are ideal for medical centres and hospitals, as they cater for complex healthcare needs. E.g. bed-bound and unconscious patients. **The Escape-Mattress** is an "A-label quality" product and produced in-house.

Our mattresses are CE certified & TÜV certified up to 150kg. We guarantee the highest quality standards.

The **Escape-Mattress Solo** is fastened to the wall. By opening the zip at the underside of the cover, the **Escape-Mattress Solo** slides out and is ready to use. Suitable for places where there is limited space, such as operating theatres.

The **Escape-Mattress** is made of fire resistant material that is resistant to blood, urine, acids and alkaline as well.

CE certified and TÜV tested and certified for 150kg payload.

There is no better and more durable quality available on the market than any of our **Escape-Mattress** range.

#### Advantages Escape Mattress Solo:

- Supplied with a storage cover
- Can be suspended from the wall
- 4 handles at each side

#### Protection of evacuee/patient:

- 2 body flaps – 1 per side
- Foot bag
- Headrest

#### Specifications Escape Mattress Solo

Article number:	21020
Free demonstration and expert advice in advance:	✓
Ready to use within a few seconds:	✓
Heat and flame resistant:	✓
Dimensions (length x width):	200 x 68 cm
Also available in a 55cm wide version:	✓
Mattress thickness:	8 cm
Including storage:	bag or case
Maximum carrying capacity:	150 kg
Space-saving:	Wall mounted
Weight:	9.5 kg incl. cover
Adjustable headrest with immobilisation band:	✓
Warranty:	2 years
4-language manual:	✓
TÜV & CE certified:	✓

#### Usage:

The **Escape-Mattress Solo** is ideal for healthcare environments like operating theatres.



Escape Mobility USA.

1261 S Lyon Street Suite 405 Santa Ana, CA 92705, USA

T +1 (855) 470-4410

W [www.escape-mobility-usa.com](http://www.escape-mobility-usa.com)

