



Escape-Sheet Economy Line



The **Escape-Sheet** is positioned permanently between the mattress and the bed frame. In the event of an emergency, the safety straps of the **Escape-Sheet** are secured over the mattress and the patient.

The patient is immobilised safely for transporting down the stairs. The **Escape-Sheet** is used frequently in hospitals and nursing homes.

The **Escape-Sheet** is widely used by hospitals and healthcare facilities. Assists in moving a patient to a place of safety, horizontally as well as vertically. Our sheets are water repellent, anti-static and also available in heat and flame resistant.

Advantages of the Escape Sheet Economy Line:

- Plastic buckle pockets increase the life expectancy of the product (allocated on the side of the sheet).
- Handgrips to assist with sideways manoeuvrability.
- Suitable for each mattress height.
- TÜV tested.
- CE certified.

Evacuation directly.... from a bed

What's an Escape Sheet?

People who aren't able to leave their bed independently need to be evacuated using evacuation equipment. One of these aids is the **Escape-Sheet**. A mattress with an **Escape-Sheet** underneath can be transformed very quickly into an evacuation mattress.

How does an Escape Sheet work?

In case of emergency the strong safety belts are taken from underneath the mattress. The belts are strapped over the patient and bedding so the patient is safely secured to the mattress. The **Escape-Sheet** can then be easily pulled to the floor and taken to a place of safety.

Our sheets are CE certified & TÜV certified up to 250kg. We guarantee the highest quality standards.

Specifications Escape Sheet Standard	
Article number:	21270
Dimensions (length x width):	± 200 x 90 cm
Material:	67% Polyester / 33% Cotton
Maximum payload:	200 kg
Weight:	ca. 800 gram
Demonstration and advise provided:	✓
Ready within seconds:	✓
Immediate evacuation directly from the bed:	✓
Industrial wash 60° C:	✓
For all mattress thicknesses:	✓
Buckles in pockets:	✓
Warranty:	2 years

Our advice:

Make sure every bed in your organization is equipped with an **Escape-Sheet**, so you're always prepared for optimal and rapid evacuation.



Escape Mobility USA.

1261 S Lyon Street Suite 405 Santa Ana, CA 92705, USA

T +1 (855) 470-4410

W www.escape-mobility-usa.com

