



# Escape-Sheet Standard



## Best choice for hospitals and healthcare facilities

### What's an Escape Sheet?

People who are unable to leave their beds, must be evacuated during an emergency. One way of evacuation is an Escape-Sheet®. An **Escape-Sheet** is permanently under the mattress and in an event of an emergency, the safety straps of the **Escape-Sheet** are secured over the mattress and the patient. The patient is then secure and immobilized for safe horizontal and vertical evacuation down the stairs.

The **Escape-Sheet** is an "A-label quality" product and produced in-house.

**Our sheets are CE certified & TÜV certified up to 250kg. We guarantee the highest quality standards.**

The **Escape-Sheet** is positioned permanently between the mattress and the bed frame. In the event of an emergency, the safety straps of the **Escape-Sheet** are secured over the mattress and the patient.

The patient is immobilised safely for transporting down the stairs. The **Escape-Sheet** is used frequently in hospitals and nursing homes.

The **Escape-Sheet** is widely used by hospitals and healthcare facilities. Assists in moving a patient to a place of safety, horizontally as well as vertically. Our sheets are water repellent, anti-static and also available in heat and flame resistant.

### Advantages of the Escape Sheet Standard:

- Straps and handles highly visible in bright green colour.
- Extra-long towing straps.
- Plastic buckle pockets increase the life expectancy of the product
- Handgrips to assist with sideways manoeuvrability.
- Static-Control™ (anti-static polyester).
- Adjusting corner straps (suitable for each mattress height).
- Hydro-tech Nano technology (water-repellent)

### Specifications Escape Sheet Standard

Article number:	21240
Dimensions (length x width):	200 x 90 cm
Material:	Polyester
Maximum payload:	250 kg
Weight:	1000 gram
Demonstration and advise provided:	✓
Ready to use within a few seconds:	✓
Immediate evacuation directly from the bed:	✓
Industrial wash 75° C:	✓
Hydro-tech nano technology (water-repellent):	✓
Adjustable corner straps:	✓
For all mattress thicknesses:	< 23 cm
Buckles in pockets:	✓
Heat and flame resistant version (HFR):	Optional
Towing straps available in 3 different colors:	Optional
Warranty:	2 Years

### Usage:

The **Escape-Sheet** is widely used by hospitals and healthcare facilities. Assists in moving a patient to a place of safety.



Escape Mobility USA.

1261 S Lyon Street Suite 405 Santa Ana, CA 92705, USA

T +1 (855) 470-4410

W www.escape-mobility-usa.com

